Tubes Radiogènes à Anode Tournante  
 Drehanoden - Röntgenröhre  
 Tubos de Rayos-X con Ánodo Giratorio

Note: Document originally drafted in the English language.

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
Manufactured by Varian Medical Systems  
Fabrique par Varian Medical Systems  
Hergestellt von Varian Medical Systems  
Fabricado por Varian Medical Systems

Specifications subject to change without notice.  
Spécifications susceptibles d'être modifiées sans préavis.  
Technische Daten ohne Gewähr.  
Especificaciones sujetas a cambio sin previo aviso.

Caractéristiques d'Émission du Filament CEI 60613  
Kathoden - Emissionskennlinien IEC 60613  
Características de Emisión del Catodo IEC 60613

**Three Phase Emission (±.15 A)**  
**GS-2076 0.9W X 0.7L** ■

Le Gaine B-205H

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Encaje de B-205H

Maximum Peak Voltage .....	140 kV
Anode to Ground .....	70 kV
Cathode to Ground .....	70 kV
Maximum X-ray Tube Assembly Heat Content .....	2.1 MJ (3.0 MHU)
Maximum Continuous Heat Dissipation (max. housing temperature 78°C) (Includes stator heat).....	3.5 kW (4.9 kHU/sec)
Maximum Heat Exchanger Dissipation .....	3.5 kW (4.9 HU/sec)
Beam Collimation* X-Ray Tube Assembly Permanent filtration .....	0.7 mm AL equivalent IEC 60522 (plus a .2 mm titanium filter attached to collimator)
1 black dot must be visible when port mounting bracket is installed	
Loading Factors for Leakage Radiation .....	140 kV, 25 mA
Weight: Tube & Housing .....	72 kg (158 lbs)
IEC Classification .....	Class 1
Overtemperature Switch - (SPST) Normally Closed ..... Opening at 68°C to 79°C, located on Oil Line Boss 5A @ 125 Vac or 2.5A @250 Vac	
Pressure Switch - Normally Closed .....	Opens at 5 PSIG ± 3 PSIG 5A @ 250 Vac
* Collimator adapter brackets in factory torqued to 15 ft. lbs. DO NOT overtorque. Overtorquing will cause compression of vibration damping material which will cause misalignment of the filter holder.	

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


**INSTALLATION DIRECTIONS**

When using a Varian GS-2076K in the B-205H Housing on a GE scanner, the installer must set the STAND-BY and BOOST currents so that the installer can successfully run "AUTOGEN-CAL". The installer should set the filaments to these values:

STAND-BY = 360 ±5 mA  
BOOST = 650 ±15 mA  
(Ref.: G.E. X-ray Service Manual (SM9803B))

The weight of the Varian housing and the GE housing are the same. The counterweights DO NOT have to be removed when installing the GS-2076K/B205H.

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VEILLE = 360 ±5 mA  
AMPLIFICATION = 650 ±15 mA  
(Ref.: G.E. X-ray Service Manual (SM9803B))

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
**TUBE HOUSING ASSEMBLY HEATING AND COOLING**

**Note:**  
1. Heat inputs into housing include tube power, filament power, and stator power.  
2. Heating curves based on no restrictions of natural convection around tube housing assembly.  
3. Heating and cooling curves reflect maximum tube performance. Tube operation is ultimately limited by system software control.

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


<b>Stator Type:</b>	
<b>Stator Coil Resistance:</b>	
Black to White	46 Ohms ±15%
Green to White	23 Ohms ±15%
Black to Green	23 Ohms ±15%
Anode rotation speed - 8000 RPM minimum	
<b>Starter Voltage:</b>	<b>Start Run</b>
150/180 Hz	385 VAC 85 VAC
<b>Time to Full Speed:</b>	
150/180 Hz	0 - 8000 RPM 10 Sec.
<b>X-Ray Tube Assembly:</b>	
GS-2076K/B-205H IEC 60601-2-28	

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**Note:**

Heating and cooling curves reflect maximum tube performance. Tube operation is ultimately limited by system software control.

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