Tubes Radiogènes à Anode Tournante  
Drehanoden - Röntgenröhre  
Tubos de Rayos-X con Ánodo Giratorio

## GS-3576S/MX165NP Housing

Note: Document originally drafted in the English language.

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
Manufactured by Varian Medical Systems  
Fabrique par Varian Medical Systems  
Hergestellt von Varian Medical Systems  
Fabricado por Varian Medical Systems


Specifications subject to change without notice.  
Spécifications susceptibles d'être modifiées sans préavis.  
Technische Daten ohne Gewähr.  
Especificaciones sujetas a cambio sin previo aviso.

Caractéristiques d'Émission du Filament CEI 60613  
Kathoden - Emissionskennlinien IEC 60613  
Características de Emisión del Catodo IEC 60613

THREE PHASE EMISSION (± .15 A)  
GS-3576S 0.9 x 0.5

THREE PHASE EMISSION (± .15 A)  
GS-3576S 1.2 x 1.2

Le Gaine MX165NP

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Encaje de MX165NP

Maximum Peak Voltage .....	140 kV
Anode to Ground .....	70 kV
Cathode to Ground .....	70 kV
Maximum X-ray Tube Assembly Heat Content .....	2.6 MJ (3.7 MHU)
Maximum Continuous Heat Dissipation (max. housing temperature 78°C) (Includes stator heat) .....	3.7 kW (5.0 kHU/seg)
Focal Point Position (Central Ray) Within 1mm (X, Y Direction from the center of radiation port.)	
X-Ray Tube Assembly	
Permanent filtration .....	1.2 mm Al IEC 60522
Loading Factors for Leakage Radiation .....	140 kV, 26 mA
Ambient Air Temperature Limits for Operation .....	5°C to 35°C
Temperature Limits for Storage and Transport .....	-34°C to +60°C
Humidity .....	10% to 90%
Atmospheric Pressure Range .....	70 kPa to 106 kPa
Weight - Tube Assembly (Tube & Housing) .....	70 kg (154 lbs)
IEC Classification .....	Class 1
Safety Devices - Thermal Switch, Normally Closed Contact .....	Opening at 74°C to 79°C
Pressure Switch Normally Closed Contact .....	Opening at 5 PSIG ±3 PSIG

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
**GS-3576S**  
**TUBE HOUSING ASSEMBLY HEATING AND COOLING**

**Note:**

- Heat inputs into housing include tube power, filament power, and stator power.
- Heating curves based on no restrictions of natural convection around tube housing assembly.
- Heating and cooling curves reflect maximum tube performance. Tube operation is ultimately limited by system software control.

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

**Housing Connection  
Stator and Safety Switches  
Plug Connector**

**Heat Exchanger Connection  
Pump and Fans 115 VAC, 60 Hz  
Plug Connector**

Stator Drive Frequency Fréquence d'entraînement du stator Statorantrieb Frequenz Frecuencia de la impulsión del estator	RPM
50 Hz	2800 - 3000
60 Hz	3400 - 3600
150 Hz	8500 - 9000
180 Hz	9500 - 10,800

<b>Stator:</b>		
<b>Stator Coil Resistance:</b>		
Black to White	23Ω ±15%	
Green to White	23Ω ±15%	
Black to Green	46Ω ±15%	
<b>Starter Voltage:</b>		
50/60 Hz	<b>Start</b> 385 VAC	<b>Run</b> 75 VAC
150/180 Hz	450 VAC	85 VAC
<b>Time to Full Speed:</b>		
50/60 Hz	0 - 2700 RPM 10 Sec.	
<b>X-Ray Tube Assembly:</b>		
GS-3576S/MX165NP IEC 60601-2-28		

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
## ANODE HEATING AND COOLING CURVES

**Note:**

1. Heating and cooling curves reflect maximum tube performance. Tube operation is ultimately limited by system software control.

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