



## ARIA Patient & Accessory Verification System

ARIA's new Patient & Accessory Verification System (PAVS) provides an additional safety measure to your radiation therapy treatment process. PAVS uses barcode technology to identify patients and verify accessories prior to treatment delivery. The system provides an efficient way for therapists to confirm that each patient is matched with the appropriate, prescribed treatment plan, and to ensure that the correct accessories are used during each treatment session.

# Match the plan to the patient

With ARIA's new Patient & Accessory Verification System (PAVS), therapists can positively match the plan displayed on the treatment console with the patient and accessories inside the treatment room. PAVS streamlines daily treatments and provides an added measure of safety to the process.

## Key features

- > Electronic patient selection from queue
- > Alert notifications for accessory changes and removal
- > Beam interlock until alerts are acknowledged
- > Verification of custom blocks, Bolus and compensators
- > Compatible with Varian's 4D Console as well as select non-Varian accelerators



---

### USA Headquarters, California

Varian Medical Systems  
Palo Alto, CA  
Tel: 650.424.5700  
800.544.4636  
Fax: 650.493.5637  
[www.varian.com](http://www.varian.com)

### Headquarters Europe, Eastern Europe, Africa, Middle & Near East

Varian Medical Systems International AG  
Zug, Switzerland  
Tel: 41.41.749.8844  
Fax: 41.41.740.3340  
[info.europe@varian.com](mailto:info.europe@varian.com)